

## 293T OS8 Cell Line

Cat. No: KC-1113

Version 18092002

I. Cell Line Information.....	1
II. Background.....	1
III. Cell Line Generation.....	2
IV. Characterization using FACS.....	2
V. Application.....	2
VI. Cell Resuscitation.....	2
VII. Cell Freezing.....	3
VIII. References.....	3

### I. Cell Line Information

<b>Catalog number</b>	KC-0113
<b>Cell line name:</b>	293T OS8 Cell Line
<b>Gene ID/Accession #:</b>	NA
<b>Host cell line</b>	293T
<b>Cell type:</b>	Human embryonic kidney
<b>Description:</b>	HEK293T cell line stable expressing membrane OKT3 scFV sequence
<b>Quantity:</b>	One vial of frozen cells (5X10 <sup>6</sup> per vial)
<b>Stability:</b>	Stable in culture over a minimum of 10 passages
<b>Application:</b>	Drug screening and biological assays
<b>Freeze medium:</b>	70% DMEM + 20% FBS + 10% DMSO
<b>Propagation medium:</b>	DMEM + 10% FBS + 50ug/ml Hygromycin
<b>Selection marker:</b>	Hygromycin
<b>Morphology:</b>	Fibroblastoid cells growing as monolayer
<b>Subculture:</b>	Split saturated culture 1:4~1:5 every 2~3 days; seed out at about 1-3 x 10 <sup>5</sup> cells/ml
<b>Incubation:</b>	37 °C with 5% CO <sub>2</sub>
<b>Doubling time:</b>	Approximately 30 hours
<b>Mycoplasma status:</b>	Negative
<b>Biosafety level:</b>	1
<b>Storage:</b>	Liquid nitrogen immediately upon receiving

### II. Background

293T-OS8 stable cell was used as an artificial antigen presenting cells (APCs), OS8 could function as a membrane anchored T cell engager that directly activates TCR in T cell-based assay.

### III. Cell Line Generation

293T OS8 cell line was generated using lentiviral vector expressing a ScFV sequence of anti-human CD3 mAb OKT3 and a C-terminal domain of mouse CD8a which consist of transmembrane and cytoplasmic domains.

### IV. Characterization using FACS




Figure: Characterization of OS8 overexpressing in 293T stable clones using FACS.

### V. Application

Hybridoma or Binders of ligand screening in T cell activation assay

### VI. Cell Resuscitation

1. Prewarm culture medium (DMEM supplemented with 10% FBS and 0.5ug/ml puromycin) in a 37°C water bath.
2. Thaw the frozen vial in a 37°C water bath for 1-2 minutes.
3. Transfer the vial into biosafety cabinet, and wipe the surface with 70% ethanol.
4. Unscrew the top of the vial and transfer the cell suspension gently into a sterile centrifuge tube containing 9.0 mL complete culture medium.
5. Spin at ~ 125 x g for 5~7 minutes at room temperature, and discard the supernatant without disturbing the pellet.
6. Resuspend cell pellet with the appropriate volume of complete medium and transfer the cell suspension into a T25 culture flask.
7. Incubate the flask at 37°C, 5% CO<sub>2</sub> incubator.
8. Split saturated culture 1:4 ~ 1:5 every 2~3 days; seed out at about 1-3 x 10<sup>5</sup> cells/ml.

## VII. Cell Freezing

1. Prepare the freezing medium (70% DMEM + 20% FBS + 10% DMSO) fresh immediately before use.
2. Keep the freezing medium on ice and label cryovials.
3. Harvest cells to a sterile, conical centrifuge tube during the logarithmic growth period and count the cells.
4. Centrifuge the cells at 250 x g for 5 minutes at room temperature and carefully aspirate off the medium.
5. Resuspend the cells at a density of at least  $3 \times 10^6$  cells/ml in chilled freezing medium.
6. Aliquot 1 ml of the cell suspension into each cryovial.
7. Freeze cells in the CoolCell freezing container overnight in a  $-80^{\circ}\text{C}$  freezer.
8. Transfer vials to liquid nitrogen for long-term storage.

## VIII. References

1. Latouche, Jean-Baptiste; Sadelain, Michel (2000). "Induction of human cytotoxic T lymphocytes by artificial antigen-presenting cells". *Nature Biotechnology*. 18 (4): 405–409.
2. Perica, Karlo; Kosmides, Alyssa K; Schneck, Jonathan P (2015). "Linking form to function: Biophysical aspects of artificial antigen presenting cell design". *Biochimica et Biophysica Acta (BBA) - Molecular Cell Research*. 1853 (4): 781–790

**Kyinno Biotechnology Co., Ltd.**

Yizhuang Biomedical Park, No. 88, Beijing, China

Tel: +86-10-58222702

E-mail: [bd@kyinno.com](mailto:bd@kyinno.com)Web: [www.kyinno.com](http://www.kyinno.com)**For Research Use Only**